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Para Leishmaniasis Cutánea

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Referencia

1. Martin SK, Thuita-Harun L, Adoyo-Adoyo M, Wasunna KM. 1998. A diagnostic ELISA for visceral leishmaniasis, based on antigen from media conditioned by *Leishmania donovani* promastigotes. *Ann Trop Med Parasitol*. 92: 571-7
2. Ryan JR, Smithyman AM, Stiteler JM, Rowton ED, Valli L, Orndorff GR, Seaman J, Thuita-Harun L and Martin SK. An ELISA based on soluble promastigote antigen detects IgM and IgG antibodies in visceral and cutaneous leishmaniasis. Abstract to November, 1998
3. Ryan JR, Rajasekariah GH, Thuita-Harun L, Yi L, Smithyman AM, Kinnamon KE, and Martin SK. The potential use of antigens shed, excreted and secreted by *Leishmania donovani* promastigotes for antigen detection tests for leishmaniasis. Abstract to December, 1999 American Society for Tropical Medicine and Hygiene Meeting, Washington, DC
4. Adoyo MA, Magiri CG, Ryan JR, Martin SK, Mason CJ and Wasunna KM. Preliminary applications of an antibody capture ELISA to monitor leishmaniasis therapy. Abstract to December, 1999 American Society for Tropical Medicine and Hygiene Meeting, Washington, DC
5. Mimori T, Matsumoto T, Calvopina MH, Gomez EA, Saya H, Katakura K, Nonaka S, Shamsuzzaman SM and Hashiguchi Y. Usefulness of sampling with a cotton swab for PCR diagnosis of cutaneous leishmaniasis in the New World. *Acta Trop* 2002; 81(3): 197-202.
6. Marques MJ, Volpini AC, Genari O, Mayrink W and Romanha AJ. Simple form of clinical sample preservation and DNA extraction from human lesions for the diagnosis of American cutaneous leishmaniasis via polymerase chain reaction. *Am J Trop Med Hyg* 2001; 65(6): 902-6.
7. Mallik MK, Das DK and Haji BE. Fine needle aspiration cytology diagnosis of cutaneous leishmaniasis in a 6-month-old child: a case report. *Acta Cytol* 2001; 45(6): 1005-7.
8. Ramirez JR, Aguledo S, Muskus C, Alzate JF, Berberich C, Barker D and Velez ID. Diagnosis of cutaneous leishmaniasis in Columbia: the sampling site within lesions influences the sensitivity of parasitologic diagnosis. *J Clin Microbiol*. 2000; 38(10): 3768-73.
9. Matsumoto T, Hashiguchi Y, Gomez EA, Calvopina MH, Nonaka S, Saya H and Mimori T. Comparison of PCR results using scrape/exudates, syringe-sucked fluid and biopsy samples

for diagnosis of cutaneous leishmaniasis in Ecuador. *Trans R Soc Trop Med Hyg.* 1999; 93(6):606-7.

10. Mohareb EW, Hanafi HA, Mikhail EM, Presley SM and Batchelor R. Evaluation of an indirect immunofluorescence assay for the detection of *Leishmania* promastigotes and amastigotes in sand flies and lesion fluid aspirates. *J Egypt Soc Parasitol.* 1998; 28(2): 313-21

11. Webb JR, Campos-Neto A, Owendale PJ, Martin TI, Stromberg EJ, Badaro R, Reed SG. Human and murine immune responses to a novel *Leishmania major* Immun. 1998; 66(7): 3279-89.

12. Rafati S, Salmanian AH, Hashemi K, Schaff C, Belli S and Fasel N. Identification of *Leishmania major* cysteine proteases as targets of the immune response in humans. *Mol Biochem Parasitol.* 2001; 113(1) 35-43.

13. de Ibarra AA, Howard JG and Snary D. Monoclonal antibodies to *Leishmania tropica major*: specificities and antigen location. *Parasitology* 1982; 85(3): 523-31.

14. Greenblatt CL, Slutzky GM, de Ibarra AA and Snary D. Monoclonal antibodies for serotyping *Leishmania* strains. *J Clin Microbiol.* 1983; 18(1): 191-3.

15. Handman E, Greenblatt CL and Goding JW. An ampipathic sulphated glycoconjugate of *Leishmania*: characterization with monoclonal antibodies. *EMBO J* 1984; 3(10): 2301-6.

16. Hailu A. The use of direct agglutination test (DAT) in serological diagnosis of Ethiopian cutaneous leishmaniasis. *Diag Microbiol and Inf Dis* 2002; 42: 251-256.

17. Beena KR, Ramesh V and Mukherjee A. Identification of parasite antigen, correlation of parasite density and inflammation in skin lesions of post kala-azar dermal leishmaniasis. *J Cutan Pathol* 2003; 30: 616-620.

18. Chargui N, Bastien P, Kallel K et al Usefulness of PCR in the diagnosis of cutaneous leishmaniasis in Tunisia. *Trans Roy Soc Trop Med and Hyg* 2005; 99, 762-768.

19. Robinson Ramirez J, Agudelo S, Muskus C et al Diagnosis of cutaneous leishmaniasis in Columbia: the sampling site within lesions influences the sensitivity of parasitologic diagnosis. *J. Clin Micro.* 2000; 38(10): 3768-3773.

20. Kassi M, Tareen I, Qazi A and Kasi PM. Fine needle aspiration cytology in the diagnosis of cutaneous leishmaniasis. *Ann Saudi Med.* 2004; 24(2): 93-97.

21. Sharquie KE, Hassen AS, Hassan SA, Al-Hamami IA. Evaluation of diagnosis of cutaneous leishmaniasis by direct smear, culture and histopathology. *Saudi Med J.* 2002; 23(8): 925-928.

22. Navin TR, Arana FE, deMerida AM, Arana BA, Castillo AL and Silvers DN. Cutaneous Leishmaniasis in Guatemala: A comparison of Diagnostic methods. *Am J Trop Med Hyg.* 1990; 42(1): 36-42.

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InBios International, Inc.